TECHNICAL SPECIFICATIONS

Weight		77-100 lbs 35-45 kg		101-121 lbs 46-55 kg		122-148 lbs 56-67 kg		149-172 lbs 68-78 kg		173-199 lbs 79-90 kg		200-242 lbs 91-110 kg	
Activity level		1/2	3	1/2	3	1/2	3	1/2	3	1/2	3	1/2	3
Audacieuse *HL01__*	> 24 cm*			XS	S	S	M	M	L	L	XL	XL	XL
Intrepide *HL02__*	> 19 cm*	PXS	PXS	PXS	PS	PS	PM	PM	PL	PL	PXL	PXL	PXL

Activity level: 1 = Low intensity / 2 = Medium intensity / 3 = Very high intensity (sprint)

***Minimum clearance**: distance between the ground and the bottom of the socket on the walking prosthesis with shoe.

FOR BILATERAL PATIENTS: Increase by one category for added stability.

Dare to dream



BORIS, TRAIL RUNNING

"This blade has the enormous advantage of being versatile. Whether you're a non-athlete, a regular or an elite athlete, you'll enjoy running with it on all types of terrain."

SARAH, HIKING AND ROADS

"Don't be afraid to challenge yourself, to have no limits. When you try, you realize that you're capable of so much!"

HOPPER

US Headquarters | 1236 West Southern Ave, Suite 101,
Tempe, AZ 85282 855.450.7300 us.proteor.com

PROTEOR

HOPPER

MORE VERSATILE THAN A TYPICAL RUNNING BLADE.
MORE DYNAMIC THAN A SPORT FOOT.



EXCLUSIVELY DISTRIBUTED BY PROTEOR

HOPPER BLADE

Engineered with the **highest quality materials** for the utmost **versatility**, **durability** and **comfort**.

THE HOPPER BLADE - elevated by AIRBUS™

- **Lightweight**, **ultra-flexible**, **dynamic carbon fiber** blade
  - Upcycled from **higher-quality** aerospace-grade material
- **Five levels of stiffness** for diverse types/levels of activity
  - **Returns energy** smoothly and comfortably
- **C-shape** allows **heel-to-toe movement** and **rolling motion**
  - Like a **natural running movement**
- **Built to last**, maintains shape over time



THE ADAPTOR

Easy to adjust with a universal top adaptor

- Both male and female included

THE OUTSOLE - enhanced by SALOMON™

- Patented sole specifically designed to **grip like a shoe**
  - Optimal for **all terrains**
- Enhances **the step through,** providing a **natural stride** feeling
- Maximizes ground contact surface, **improving stability**
- **Interchangeable,** easy to remove/clean/change

Built for versatile use, on flat or uneven terrain

Used for:
• Running • Dancing • Hiking • Badminton • Tennis • Many other sports



THE FOUNDER'S WORDS

"My dream was to be able to run with my children, but I was confronted with **different issues.**
This raised the question: How to make sport accessible to amputees?"

Jérôme,
co-founder of Hopper

TWO MODELS, DIFFERENT CLEARANCES

AUDACIEUSE

Audacious

Since 2021, the pioneering **Audacieuse** has been helping people around the world rediscover the exhilarating joy of sporting activities.

Minimum clearance: 24 cm.



INTREPIDE

Intrepid

The **Intrepide**, recently developed upon request from users (with feedback from prosthetists) is designed for adventurers of all ages who want to push limits with boundless energy.

Minimum clearance: 19 cm.

OUR MISSION

Sport isn't reserved for elite athletes, nor is it just a casual hobby. It is essential for health, well-being, and social inclusion. **Our mission is to make sport accessible to everyone, regardless of performance or financial considerations.**